


Supporting Information

Fig S1. Alignments of sequence records derived from two type specimens of Geometridae, one with high quality DNA (a) and one with low quality DNA (b). The alignments show only a single representative of each distinct sequence. In many cases, there were hundreds or thousands of a particular sequence. High quality reads have high coverage across the entire 658bp barcode region and originate from a single source – indicated by a single nucleotide (color) at each position in the contig. Low quality reads do not span the entire barcode region (i.e. they have regions lacking coverage) and often originate from multiple sources – indicated by multiple nucleotides (colors) at certain positions in the contig.

10 *See associated image file “Figure S1.svg”*



18 S1 Table. Type specimens analyzed, including sequencing results and accession numbers.

Process ID (Sanger/NGS)	Age (Yrs)	Identification	Status	Sanger Group	No. NGS Reads	Min. Cov.	Max. Cov.	Avg. Cov.	Recovered bp by NGS	NGS Contig GenBank Acc.	Sequence Read Archive Acc.
PNGTY381-13 / PNGTY1837-15	104	<i>Myrioblephara mixticolor</i>	Syntype	HQ	143804	7	115751	29924	658	KR070780	SRR1867808
PNGTY404-13 / PNGTY1827-15	112	<i>Cassephyra plenimargo</i>	Holotype	HQ	213007	72	146189	42992	658	KR070771	SRR1867811
PNGTY417-13 / PNGTY1843-15	109	<i>Psilalcis auropurpurea</i>	Syntype	HQ	106286	0	74012	20477	448	KR070767	SRR1867812
PNGTY466-13 / PNGTY1839-15	110	<i>Paralcidia marginata</i>	Syntype	HQ	221885	5	168474	44541	657	KR070779	SRR1867813
PNGTY473-13 / PNGTY1823-15	110	<i>Atmoceras plumosa</i>	Syntype	HQ	143340	30	76376	28215	658	KR070768	SRR1867814
PNGTY1047-14 / PNGTY1845-15	110	<i>Tripteridia viridisecta</i>	Syntype	HQ	188107	1	101191	37855	570	KR070783	SRR1867815
PNGTY1070-14 / PNGTY1834-15	111	<i>Gymnoscelis ochriplaga</i>	Holotype	HQ	186897	1	103399	38169	657	KR070772	SRR1867816
PNGTY1098-14 / PNGTY1824-15	110	<i>Axinoptera fasciata</i>	Holotype	HQ	166116	0	83389	31838	474	KR070785	SRR1867817
PNGTY1106-14 / PNGTY1826-15	112	<i>Calluga semirasata</i>	Holotype	HQ	215946	1	154302	43408	658	KR070769	SRR1867818
PNGTY1124-14 / PNGTY1833-15	123	<i>Eois semirubra</i>	Holotype	HQ	232024	106	143908	46803	658	KR070773	SRR1867819
PNGTY142-12 / PNGTY1831-15	112	<i>Collix ghosha dichobathra</i>	N/A	MQ	11665	0	11082	2142	459	KR070778	SRR1945335
PNGTY158-12 / PNGTY1838-15	111	<i>Papuarisme brunneata</i>	Holotype	MQ	6479	0	5747	1165	569	KR070787	SRR1945382
PNGTY189-12 / PNGTY1835-15	109	<i>Hyposidra apicifulva</i>	N/A	MQ	62208	0	45441	12301	570	KR070781	SRR1945383
PNGTY801-13 / PNGTY1836-15	101	<i>Milionia knowlei</i>	Syntype	MQ	44190	1	43301	9153	658	KR070777	SRR1945384
PNGTY587-13 / PNGTY1846-15	118	<i>Ctimene basistraga obsoleta</i>	Syntype	MQ	546	0	105	31	323	KR070784	SRR1946575
PNGTY639-13 / PNGTY1842-15	102	<i>Pseudeusemia bursadoides dignitosa</i>	Syntype	MQ	134542	37	82031	27382	658	KR070774	SRR1945385
PNGTY917-13 / PNGTY1840-15	108	<i>Pingasa nobilis furvifrons</i>	Holotype	MQ	46793	6	24276	9516	658	KR070782	SRR1945386
PNGTY923-13 / PNGTY1821-15	121	<i>Aeolochroma caesia</i>	Holotype	MQ	68837	3	43442	14002	658	KR070764	SRR1945387
PNGTY957-13 / PNGTY1844-15	106	<i>Sarcinodes subvirgata</i>	Holotype	MQ	99655	86	42405	20379	658	KR070762	SRR1945388
PNGTY971-13 / PNGTY1828-15	105	<i>Celerena lerne amplimargo</i>	Holotype	MQ	113363	1	42333	21657	569	KR070763	SRR1945389
PNGTY475-13 / PNGTY1832-15	104	<i>Dyscheralcis retroflexa</i>	Syntype	LQ	2681	0	1278	424	514	KR070766	SRR1867935
PNGTY1146-14 / PNGTY1822-15	110	<i>Alcis irrufata</i>	Holotype	LQ	49944	2	37401	8881	657	KR070770	SRR1867936

PNGTY1155-14 / PNGTY1830-15	120	<i>Cleora repetita</i> ab. <i>suffusa</i>	Holotype	LQ	7632	1	5708	731	454	KR070765	SRR1867937
PNGTY008-12* / N/A	109	<i>Spectrobasis</i> <i>differens</i>	Syntype	LQ	1468	0	1157	280	237	N/A	SRR1867938
PNGTY073-12* / N/A	104	<i>Desmoclystia</i> <i>unipuncta</i>	Syntype	LQ	320	0	116	40	357	N/A	SRR1867939
PNGTY102-12* / N/A	120	<i>Sterrhochaeta</i> <i>minuta</i>	Syntype	LQ	3081	0	215	91	324	N/A	SRR1867940
PNGTY120-12* / N/A	118	<i>Propitex alternata</i>	Holotype	LQ	14863	0	1554	163	323	N/A	SRR1867941
PNGTY756-13 / PNGTY1825-15	105	<i>Bursadopsis</i> <i>plenifascia</i>	Syntype	LQ	133263	4	71444	20742	634	KR070786	SRR1867942
PNGTY1080-14 / PNGTY1829-15	117	<i>Chloroclystis</i> <i>rufofasciata</i>	Syntype	LQ	4411	0	1685	647	419	KR070775	SRR1867943
PNGTY1128-14 / PNGTY1841-15	112	<i>Polyacme</i> <i>straminea</i> ab. <i>brunneata</i>	Holotype	LQ	141402	23	71197	28117	658	KR070776	SRR1867944

19 The four Process ID's marked with an asterisk (*) represent specimens where NGS analysis generated
20 sequence reads from multiple species. HQ – high quality; MQ – medium quality; LQ – low quality; N/A –
21 not applicable.
22
23
24 **Table S2. Primers used in the first (PCR1) and second (PCR2) reactions to allow the analysis of 10**
25 **specimens in an Ion Torrent PGM run.**
26

PCR	Code	Primer Name	Sequence (5'-3')	MID	Adapter
PCR1	F1	ARTH-NGS-F1.1-ion1	CTAAGGTAACATTCAACCAATCATAAGATATTGG	None	None
	F2	ARTH-NGS-F2.1-ion1	CTAAGGTAACATRRWRATGATCAARTWTATAAT	None	None
	F3	ARTH-NGS-F3.1-ion1	CTAAGGTAACTTATAATTGGDGGRTTTGGWAATTG	None	None
	F4	ARTH-NGS-F4.1-ion1	CTAAGGTAACAGWAGWATWRTWRAWAVWGG	None	None
	F5	ARTH-NGS-F5.1-ion1	CTAAGGTAACATTTTWSWCTWCATWTDGCWGG	None	None
	F6	ARTH-NGS-F6.1-ion1	CTAAGGTAACATTTGTWTGAKCWRTWKKWATTAC	None	None
PCR1	R1	ARTH-NGS-R1.1-ion1	CTAAGGTAACWGGTATWACTATRAAAAAATTAT	None	None
	R2	ARTH-NGS-R2.1-ion2	TAAGGAGAATCCTAAWCTWATRTTRTTADWCG	None	None
	R3	ARTH-NGS-R3.1-ion3	AAGAGGATTCARDGGDGGRTAWACWGTTCAWCC	None	None
	R4	ARTH-NGS-R4.1-ion4	TACCAAGATCGTWGWAATRAATTTDATWGCWCC	None	None
	R5	ARTH-NGS-R5.1-ion5	CAGAAGGAACGTTARWARTATDGTRATDGCWCC	None	None

	R6	ARTH-NGS-R6.1-ion6	CTGCAAGTCTAAACTTCTGGATGTCCAAAAATCA	None	None
PCR2	F1	ARTH-NGS-F1.2-ion1	CCATCTCATCCCTGCGTGTCTCCGACTCAGCTAAGGTAACATTCAACCAATCATAAAGATATTGG	lonXpress1	A
	F2	ARTH-NGS-F2.2-ion1	CCATCTCATCCCTGCGTGTCTCCGACTCAGCTAAGGTAACATRRWRATGATCAARTWATAAT	lonXpress1	A
	F3	ARTH-NGS-F3.2-ion1	CCATCTCATCCCTGCGTGTCTCCGACTCAGCTAAGGTAACTTATAATTGGDGGRTTTGGWAATTG	lonXpress1	A
	F4	ARTH-NGS-F4.2-ion1	CCATCTCATCCCTGCGTGTCTCCGACTCAGCTAAGGTAACAGWAGWATWRTWRAWAVWGG	lonXpress1	A
	F5	ARTH-NGS-F5.2-ion1	CCATCTCATCCCTGCGTGTCTCCGACTCAGCTAAGGTAACATTTTWSWCTWCATWTDGCWGG	lonXpress1	A
	F6	ARTH-NGS-F6.2-ion1	CCATCTCATCCCTGCGTGTCTCCGACTCAGCTAAGGTAACATTGTGTGAKCWRTWKKWATTAC	lonXpress1	A
PCR2	F1	ARTH-NGS-F1.2-ion2	CCATCTCATCCCTGCGTGTCTCCGACTCAGTAAGGAGAACATTCAACCAATCATAAAGATATTGG	lonXpress2	A
	F2	ARTH-NGS-F2.2-ion2	CCATCTCATCCCTGCGTGTCTCCGACTCAGTAAGGAGAACATRRWRATGATCAARTWATAAT	lonXpress2	A
	F3	ARTH-NGS-F3.2-ion2	CCATCTCATCCCTGCGTGTCTCCGACTCAGTAAGGAGAACTTATAATTGGDGGRTTTGGWAATTG	lonXpress2	A
	F4	ARTH-NGS-F4.2-ion2	CCATCTCATCCCTGCGTGTCTCCGACTCAGTAAGGAGAACAGWAGWATWRTWRAWAVWGG	lonXpress2	A
	F5	ARTH-NGS-F5.2-ion2	CCATCTCATCCCTGCGTGTCTCCGACTCAGTAAGGAGAACATTTTWSWCTWCATWTDGCWGG	lonXpress2	A
	F6	ARTH-NGS-F6.2-ion2	CCATCTCATCCCTGCGTGTCTCCGACTCAGTAAGGAGAACATTGTGTGAKCWRTWKKWATTAC	lonXpress2	A
PCR2	F1	ARTH-NGS-F1.2-ion3	CCATCTCATCCCTGCGTGTCTCCGACTCAGAAGAGGATTCATTCAACCAATCATAAAGATATTGG	lonXpress3	A
	F2	ARTH-NGS-F2.2-ion3	CCATCTCATCCCTGCGTGTCTCCGACTCAGAAGAGGATTCATRRWRATGATCAARTWATAAT	lonXpress3	A
	F3	ARTH-NGS-F3.2-ion3	CCATCTCATCCCTGCGTGTCTCCGACTCAGAAGAGGATTCATTATAATTGGDGGRTTTGGWAATTG	lonXpress3	A
	F4	ARTH-NGS-F4.2-ion3	CCATCTCATCCCTGCGTGTCTCCGACTCAGAAGAGGATTCAGWAGWATWRTWRAWAVWGG	lonXpress3	A
	F5	ARTH-NGS-F5.2-ion3	CCATCTCATCCCTGCGTGTCTCCGACTCAGAAGAGGATTCATTTTWSWCTWCATWTDGCWGG	lonXpress3	A
	F6	ARTH-NGS-F6.2-ion3	CCATCTCATCCCTGCGTGTCTCCGACTCAGAAGAGGATTCATTGTGTGAKCWRTWKKWATTAC	lonXpress3	A
PCR2	F1	ARTH-NGS-F1.2-ion4	CCATCTCATCCCTGCGTGTCTCCGACTCAGTACCAAGATCATTCAACCAATCATAAAGATATTGG	lonXpress4	A
	F2	ARTH-NGS-F2.2-ion4	CCATCTCATCCCTGCGTGTCTCCGACTCAGTACCAAGATCATRRWRATGATCAARTWATAAT	lonXpress4	A
	F3	ARTH-NGS-F3.2-ion4	CCATCTCATCCCTGCGTGTCTCCGACTCAGTACCAAGATCTTATAATTGGDGGRTTTGGWAATTG	lonXpress4	A
	F4	ARTH-NGS-F4.2-ion4	CCATCTCATCCCTGCGTGTCTCCGACTCAGTACCAAGATCAGWAGWATWRTWRAWAVWGG	lonXpress4	A
	F5	ARTH-NGS-F5.2-ion4	CCATCTCATCCCTGCGTGTCTCCGACTCAGTACCAAGATCATTTTWSWCTWCATWTDGCWGG	lonXpress4	A
	F6	ARTH-NGS-F6.2-ion4	CCATCTCATCCCTGCGTGTCTCCGACTCAGTACCAAGATCTATTGTGTGAKCWRTWKKWATTAC	lonXpress4	A
PCR2	F1	ARTH-NGS-F1.2-ion5	CCATCTCATCCCTGCGTGTCTCCGACTCAGCAGAAGGAACATTCAACCAATCATAAAGATATTGG	lonXpress5	A
	F2	ARTH-NGS-F2.2-ion5	CCATCTCATCCCTGCGTGTCTCCGACTCAGCAGAAGGAACATRRWRATGATCAARTWATAAT	lonXpress5	A
	F3	ARTH-NGS-F3.2-ion5	CCATCTCATCCCTGCGTGTCTCCGACTCAGCAGAAGGAACTTATAATTGGDGGRTTTGGWAATTG	lonXpress5	A
	F4	ARTH-NGS-F4.2-ion5	CCATCTCATCCCTGCGTGTCTCCGACTCAGCAGAAGGAACAGWAGWATWRTWRAWAVWGG	lonXpress5	A
	F5	ARTH-NGS-F5.2-ion5	CCATCTCATCCCTGCGTGTCTCCGACTCAGCAGAAGGAACATTTTWSWCTWCATWTDGCWGG	lonXpress5	A
	F6	ARTH-NGS-F6.2-ion5	CCATCTCATCCCTGCGTGTCTCCGACTCAGCAGAAGGAACATTGTGTGAKCWRTWKKWATTAC	lonXpress5	A
PCR2	F1	ARTH-NGS-F1.2-ion6	CCATCTCATCCCTGCGTGTCTCCGACTCAGCTGCAAGTTCATTCAACCAATCATAAAGATATTGG	lonXpress6	A

	F2	ARTH-NGS-F2.2-ion6	CCATCTCATCCCTGCGTGTCTCCGACTCAGCTGCAAGTTCATRRWRATGATCAARTWTA TAAT	IonXpress6	A
	F3	ARTH-NGS-F3.2-ion6	CCATCTCATCCCTGCGTGTCTCCGACTCAGCTGCAAGTTCATTATAATTGGDGGRTTTGGW AATTG	IonXpress6	A
	F4	ARTH-NGS-F4.2-ion6	CCATCTCATCCCTGCGTGTCTCCGACTCAGCTGCAAGTTCAGWAGWATWRTWRAWAVW GG	IonXpress6	A
	F5	ARTH-NGS-F5.2-ion6	CCATCTCATCCCTGCGTGTCTCCGACTCAGCTGCAAGTTCATTTTWSWCTWCATWTDGC WGG	IonXpress6	A
	F6	ARTH-NGS-F6.2-ion6	CCATCTCATCCCTGCGTGTCTCCGACTCAGCTGCAAGTTCATTGTGTGAKCWRTWKK WATTAC	IonXpress6	A
PCR2	F1	ARTH-NGS-F1.2-ion7	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCGTGATTCATTCAACCAATCATAAGATA TTGG	IonXpress7	A
	F2	ARTH-NGS-F2.2-ion7	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCGTGATTCATRRWRATGATCAARTWTA TAAT	IonXpress7	A
	F3	ARTH-NGS-F3.2-ion7	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCGTGATTCATTATAATTGGDGGRTTTGGW AATTG	IonXpress7	A
	F4	ARTH-NGS-F4.2-ion7	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCGTGATTCAGWAGWATWRTWRAWAVW GG	IonXpress7	A
	F5	ARTH-NGS-F5.2-ion7	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCGTGATTCATTTTWSWCTWCATWTDGC WGG	IonXpress7	A
	F6	ARTH-NGS-F6.2-ion7	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCGTGATTCATTGTGTGAKCWRTWKKW ATTAC	IonXpress7	A
PCR2	F1	ARTH-NGS-F1.2-ion8	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCGGATAACATTCAACCAATCATAAGAT ATTGG	IonXpress8	A
	F2	ARTH-NGS-F2.2-ion8	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCGGATAACATRRWRATGATCAARTWTA TAAT	IonXpress8	A
	F3	ARTH-NGS-F3.2-ion8	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCGGATAACATTATAATTGGDGGRTTTGGW AATTG	IonXpress8	A
	F4	ARTH-NGS-F4.2-ion8	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCGGATAACAGWAGWATWRTWRAWAVW GG	IonXpress8	A
	F5	ARTH-NGS-F5.2-ion8	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCGGATAACATTTTWSWCTWCATWTDGC WGG	IonXpress8	A
	F6	ARTH-NGS-F6.2-ion8	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCGGATAACTATTGTGTGAKCWRTWKK WATTAC	IonXpress8	A
PCR2	F1	ARTH-NGS-F1.2-ion9	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGAGCGGAACATTCAACCAATCATAAGAT ATTGG	IonXpress9	A
	F2	ARTH-NGS-F2.2-ion9	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGAGCGGAACATRRWRATGATCAARTWT ATAAT	IonXpress9	A
	F3	ARTH-NGS-F3.2-ion9	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGAGCGGAACATTATAATTGGDGGRTTTGGW AATTG	IonXpress9	A
	F4	ARTH-NGS-F4.2-ion9	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGAGCGGAACAGWAGWATWRTWRAWAVW GG	IonXpress9	A
	F5	ARTH-NGS-F5.2-ion9	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGAGCGGAACATTTTWSWCTWCATWTDG CWGG	IonXpress9	A
	F6	ARTH-NGS-F6.2-ion9	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGAGCGGAACATTGTGTGAKCWRTWKK WATTAC	IonXpress9	A
PCR2	F1	ARTH-NGS-F1.2-ion10	CCATCTCATCCCTGCGTGTCTCCGACTCAGCTGACCGAACATTCAACCAATCATAAGAT ATTGG	IonXpress10	A
	F2	ARTH-NGS-F2.2-ion10	CCATCTCATCCCTGCGTGTCTCCGACTCAGCTGACCGAACATRRWRATGATCAARTWTA TAAT	IonXpress10	A
	F3	ARTH-NGS-F3.2-ion10	CCATCTCATCCCTGCGTGTCTCCGACTCAGCTGACCGAACATTATAATTGGDGGRTTTGGW AATTG	IonXpress10	A
	F4	ARTH-NGS-F4.2-ion10	CCATCTCATCCCTGCGTGTCTCCGACTCAGCTGACCGAACAGWAGWATWRTWRAWAVW GG	IonXpress10	A
	F5	ARTH-NGS-F5.2-ion10	CCATCTCATCCCTGCGTGTCTCCGACTCAGCTGACCGAACATTTTWSWCTWCATWTDG CWGG	IonXpress10	A
	F6	ARTH-NGS-F6.2-ion10	CCATCTCATCCCTGCGTGTCTCCGACTCAGCTGACCGAACATTGTGTGAKCWRTWKK WATTAC	IonXpress10	A
PCR2	R1	ARTH-NGS-R1.2-ion1-trP1	CCTCTCTATGGGCAGTCGGTGATTAAGGTAACWGGTATWACTATRAAAAAATTAT	IonXpress1	trP1
	R2	ARTH-NGS-R2.2-ion2-trP1	CCTCTCTATGGGCAGTCGGTGATTAAGGAGAATCRAAWCTWATRTTTADWCG	IonXpress2	trP1
	R3	ARTH-NGS-R3.2-ion3-trP1	CCTCTCTATGGGCAGTCGGTGATAAGAGGATTCARDGGDGGRTAWACWGTTCAWCC	IonXpress3	trP1

R4	ARTH-NGS-R4.2-ion4-trP1	CCTCTCTATGGGCAGTCGGTGATACCAAGATCGTWGWAATRAARTTDDATWGCWCC	IonXpress4	trP1
R5	ARTH-NGS-R5.2-ion5-trP1	CCTCTCTATGGGCAGTCGGTGATCAGAAGGAACGTTARWARTATDGTRATDGCWCC	IonXpress5	trP1
R6	ARTH-NGS-R6.2-ion6-trP1	CCTCTCTATGGGCAGTCGGTGATCTGCAAGTCTAAACTTCTGGATGTCCAAAAATCA	IonXpress6	trP1

The “Code” column refers to primer labels in Fig. 1. The COI binding region within each primer sequence is shown in black, while the 10bp tail (PCR1) or MID tag (PCR2) is shown in blue. The “key sequence” (required for Ion Torrent sequencing) is shown in green and the sequencing adapters are shown in red. The 10bp tails on the PCR1 primers are technically IonXpress MID tags, but they serve only to block short amplicons from acting as primers during PCR1. They were chosen over random decamer tails to maximize primer-template matching in PCR2. The same forward and reverse PCR1 primers are used for all ten samples in the first round of PCR. In the second round of PCR, samples are discriminated by using ten different sets of MID-tagged forward PCR2 primers (the same set of PCR2 reverse primers is used for all ten samples).

Table S3. Components of PCR reactions in the NGS protocol.

	PCR 1.1	PCR 1.2	PCR 2.1, 2.2, 2.3, 2.4	PCR 2.5	PCR 2.6
10% Trehalose	5.125 µL	5.25 µL	5.75 µL	5.875 µL	6.0 µL
H₂O	0.13 µL	0.13 µL	0.13 µL	0.13 µL	0.13 µL
5X Buffer	2.5 µL	2.5 µL	2.5 µL	2.5 µL	2.5 µL
25 mM MgCl₂	1.25 µL	1.25 µL	1.25 µL	1.25 µL	1.25 µL
10 µM primers	0.125 µL each	0.125 µL each	0.125 µL each	0.125 µL each	0.125 µL each
10 µM dNTP	0.0625 µL	0.0625 µL	0.0625 µL	0.0625 µL	0.0625 µL
Taq (5U/ µL)	0.06 µL	0.06 µL	0.06 µL	0.06 µL	0.06 µL
Template	2 µL	2 µL	2 µL	2 µL	2 µL
TOTAL	12.5 µL	12.5 µL	12.5 µL	12.5 µL	12.5 µL

Reactions differ only in the number of primers and the amount of trehalose. Trehalose sourced from Fluka Analytical; Hyclone ultra-pure water from Thermo Fisher Scientific; Buffer (- MgCl₂), MgCl₂, and Taq polymerase from KAPA Biosystems (while standard CCDB protocols utilize Platinum Taq, KAPA Taq was found to be less prone to co-amplifying trace amounts of residual DNA - derived from the Taq manufacturing process - and is therefore more amenable to high cycle PCR); primers from Integrated DNA Technologies.